



AseptoPowder Hot

Powerful alkaline basic cleaner (AseptoPowder) plus eight sachets of chlorinated boosted cleaner (AseptoHot) for milking equipment and milk tanks.

- ▲ Chlorinated cleaning booster
- ▲ Readily soluble powder

APPLICATION

Daily application: Pre-clean with water. Clean with **AseptoPowder** at 0.5% (minimum 10 min. at 15°C - 50°C). Rinse with water of drinking quality.

Weekly application: Pre-clean with water. Clean with **AseptoPowder** and **AseptoHot** at 0.25% (minimum 10 min. at 60°C - 80°C). Rinse with water of drinking quality. To be used in daily rotation with **HorolithPowder**.



FEATURES AND BENEFITS

AseptoPowder Hot is a package of two products in one. It consists of **AseptoPowder** and **AseptoHot**. **AseptoPowder** is an alkaline cleaner for use in all kind of milking equipment, including cooling tanks. It has very good cleaning properties and is free of chlorine. **AseptoHot** is an alkaline chlorinated cleaning booster, which is used periodically, e.g. once per week, together with the daily used product, depending on individual demands of a plant. Both are readily soluble powders.

The **Powder** products are suitable for use with all water hardnesses.

SPECIFICATIONS

AseptoPowder Hot	
Physical Form	AseptoPowder: Light yellow powder AseptoHot: Light yellow powder
Active Substance	AseptoPowder: Sodium hydroxide, complexing agents AseptoHot: Sodium dichloros triazinetrione
pH at 20°C:	AseptoPowder: 11.5 - 12.5 AseptoHot: 11.5 - 12.5
Density at 20°C (g/cm ³):	AseptoPowder: 1.6 - 1.7 AseptoHot: 1.6 - 1.7
Storage temperature:	AseptoPowder: -10°C to +40°C AseptoHot: -10°C to +30°C

LOGISTICS

Packaging	11 kg
Product code	2331640
EAN code	4028162331646

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtstrasse 7
CH-8304 Wallisellen
www.ecolab.com