



HorolithCool

Powerful acidic cleaner for cold or hot cleaning of milking equipment and milk tanks.

- ▲ Hot or cold cleaning
- ▲ Hard or soft water

APPLICATION

Pre-clean with water. Clean with **HorolithCool** 1.0% (at room temperature, 10 min.) or 0.7% (at 70°C; 10 min.). Rinse with water of drinking quality. To be used once or twice a week in rotation with **AseptoCool**.



FEATURES AND BENEFITS

HorolithCool is a blend of acids which assures very good scale removing and cleaning properties. It is suitable for all kind of standard milking equipment, including cooling tanks. The **Cool** products can be used at a wide temperature range including cold washes and in all water hardnesses.

SPECIFICATIONS

HorolithCool	
Physical Form	Clear, colourless liquid
Active Substance	Inorganic acids
pH at 20°C:	1.0
Density at 20°C (g/cm ³):	1.20 - 1.22
Storage temperature:	-30°C to +40°C

LOGISTICS

Packaging	12 kg	24 kg	250 kg
Product code	2332040	2332060	2332070
EAN code	4028162332049	4028162332063	4028162332070

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtistrasse 7
CH-8304 Wallisellen
www.ecolab.com