



Exelerate™ Enzymatic CIP 100

Eliminate soils. Eliminate compromise.

Exelerate™ Enzymatic CIP 100 is a modern, targeted CIP cleaner that can help dairy processors meet both efficiency and food safety and quality targets. By eliminating the tradeoffs of alkaline alternatives, Exelerate™ Enzymatic unlocks a smarter approach to CIP.

Activate Precision Cleaning

Enzymatic chemistry targets and breaks down the soils specific to dairy operations, delivering a thorough clean that doesn't damage core equipment.

Elevate Safety

Exelerate™ Enzymatic is not hazardous at use solution. A lower pH further lowers the risk to employee safety relative to traditional caustic.

Optimize Operations

Exelerate™ Enzymatic enables faster rinse cycles that require less water.

Enhance Effluent Quality

Lower sodium levels support the mitigation of effluent challenges ranging from pH balancing to discharge limits.



Exelerate™ Enzymatic
users report
up to 25%
in rinse time savings
per CIP line.*



Up to 70%
reduction in sodium and
lower pH compared to
caustic alternatives.*

*Based on results from customer trials, which are specific to those individual customers. Results will vary for other customers based on factors and circumstances in their operations.

The end of alkaline-as-usual

Compared to alkaline alternatives, like caustic, Exelerate™ Enzymatic drives more efficient CIP, eliminates damage to core equipment, and simplifies effluent compliance — all while delivering an targeted clean.



**Formulated
for reuse**



**Easy dosing
monitoring**



**Fast rinses and
deactivation**



**Traceable and
measurable**

IDEAL ACROSS A NUMBER OF KEY APPLICATIONS



Trucks



Tanks



Lines



Fillers



Vats

SUCCESS STORY

200% Annual ROI

Annual value delivered for a US cream-cheese plant that adopted our enzymatic CIP cleaner



Safety

~\$99K serious injuries avoided



Water

625K gallons rinse
water savings
\$4.5K saved



Effluent Quality

\$40K+ value in
improved treatment



Concentrated Product Optimization (Chemistry)

7K less gallons



Productivity

167 hrs less CIP time



Asset Protection

\$35K gasket savings

Ready to think beyond caustic?

Contact your Ecolab representative today to learn more.

ecolab.com/exelerateenzymatic